

Quick Reference Guide


ThinPrep® Pap Test Instructions for specimen collection

ThinPrep Pap Test Rovers® Cervex-Brush® Combi device protocol

- Prepare all equipment before starting the procedure.
- Note expiry date on sample collection vial. Do not use expired vials.
- Ensure the entire plastic seal is removed from the lid of the vial and discarded.
- Complete patient details on both the request form and the vial. Specimens may be returned if details are missing from the vial.
- Remove the lid from the vial before taking the sample.
- If lubrication of the speculum is required, a little warm water should be used. If necessary a small amount of water-soluble carbomer-free lubricant should be applied sparingly to the outer portion of the speculum with great care to avoid the tip as it may interfere with the ThinPrep Pap Test process.



Record patient ID on the vial as required by local regulations.
The patient information and medical history on the cytology requisition form.



Obtain an adequate sample from the cervix using a Rovers Cervex-Brush Combi device (broom-like device with an integrated endocervical sampler). If desired, use lukewarm water to warm and lubricate the speculum. Sparingly apply water-soluble carbomer-free gel lubricant to the posterior blade of the speculum if necessary.¹² Insert the endocervical sampler part of the device into the endocervical canal deep enough to allow the other bristles to fully contact the ectocervix. Push gently, and rotate the brush 360° in a clockwise direction 2 times.



Rinse the Rovers Cervex-Brush Combi device immediately into the PreservCyt® Solution vial by pushing it into the bottom of the vial 10 times, forcing the bristles apart. As a final step, swirl the brush **vigorously** to further release material. Visually inspect the Rovers Cervex-Brush Combi device to ensure that no material remains attached. Discard the collection device. **Do not leave the head of the Rovers Cervex-Brush Combi device in the vial.**



Tighten the cap so that the torque line on the cap passes the torque line on the vial.
Do not over-tighten.



Place the vial and requisition in a specimen bag for transport to the laboratory.

Frequently Asked Questions

What do I do if the vial is knocked over?

If the vial is knocked over after the sample has been taken, look to see if any fluid still remains in the vial. If about one third of the fluid still remains in the vial, then the vial can still be sent to the lab for processing. Make a note on the request form to inform the laboratory of the sample spillage. Do not 'top up' the vial from a new vial as the lab will make any necessary level adjustments.

Can I take a ThinPrep sample if the patient is menstruating?

No. Peri-menstrual samples should be avoided.

Can I use the Rovers Cervex-Brush Combi device where there is stenosis of the os?

The Rovers Cervex-Brush Combi device has an integrated endocervical brush that can be used in the majority of cases. When the stenosis is severe, the use of only an endocervical brush may be considered.

Can I use the Rovers Cervex-Brush Combi device for vault samples?

No, the standard Rovers Cervex-Brush is advised for taking vault samples. The device should be applied by using a sweeping action.

Can I use the Rovers Cervex-Brush Combi device if there is an extensive ectropion area?

Yes, the Rovers Cervex-Brush Combi device can still be used as normal. But for the ectropion area an additional standard Rovers Cervex-Brush should be used by applying a circular sweeping action. Both sampling devices should be rinsed in the same vial.

Can I use a lubricant?

If lubrication of the speculum is required, a little warm water should be used. If necessary a small amount of water-soluble carbomer-free lubricant should be applied sparingly to the outer portion of the speculum with great care to avoid the tip as it may interfere with the ThinPrep Pap Test process.

Do I need to take a separate sample for HPV testing?

No. HPV testing can be performed on the same material collected for the ThinPrep Pap Test.

Do's

- Do gently remove any excessive mucus before sampling
- Do rotate the Rovers Cervex-Brush Combi device 2 times in a clockwise direction through 360°
- Do insert the endocervical sampler part of the device into the endocervical canal deep enough to allow the other bristles to fully contact the ectocervix
- Do insert the Rovers Cervex-Brush Combi device immediately into the vial by pushing it into the bottom of the vial 10 times, forcing the bristles apart. As a final step, swirl the brush vigorously to further release material. Visually inspect the Rovers Cervex-Brush Combi device to ensure that no material remains attached. Discard the collection device.
- Do replace the lid and tighten so that the black torque line on the cap passes the black torque line on the vial to avoid leakage
- Do keep the unlabelled portion of the sample vial free of labels so that the contents can be seen
- Do send samples to the laboratory without delay

Don'ts

- Do not routinely clean the cervix or take a cervical swab before taking a cervical sample
- Do not use non water-soluble lubricants containing carbomers
- Do not leave the collection device sitting in the vial whilst dealing with the patient
- Do not over tighten the lid onto the vial
- Do not detach the head of the Rovers Cervex-Brush Combi device into the vial
- Do not place multiple labels on the outside of the vial
- Do not apply any physician bar-coded labels vertically on the vial, apply them horizontally
- Do not apply any labels over the lid of the container
- Do not use expired vials
- Do not delay the sending of the vials to the laboratory
- Do not use during pregnancy

1. Cervicovaginal Cytology Based on the Papanicolaou Technique; Approved Guideline – Third Edition (Clinical and Laboratory Standards Institute GP15-A3).
2. Hologic internal study, data on file.
3. Rovers Cervex-Brush Combi, Information For Use, 2010-11.

Diagnostic Solutions | Hologic.com | euinfo@hologic.com

SS-05867-EUR-EN Rev. 002 © 2016 Hologic, Inc. All rights reserved. Specifications are subject to change without prior notice. Hologic, Science of Sure, PreservCyt, ThinPrep and associated logos are trademarks and/or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks, and product names are the property of their respective owners. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your Hologic representative or write to euinfo@hologic.com

This is a general information tool for medical professionals and is not a complete representation of the product(s)' Instruction for Use (IFU) or Package Insert, and it is the medical professionals' responsibility to read and follow the IFU or Package Insert. The information provided may suggest a particular technique or protocol however it is the sole responsibility of the medical professional to determine which technique or protocol is appropriate. At all times, clinicians remain responsible for utilizing sound patient evaluation and selection practices, and for complying with all applicable rules and regulations regarding accreditation, anesthesia, reimbursement, and all other aspects of in-office procedures. In no event shall Hologic be liable for damages of any kind resulting from your use of the information presented.