



### **Non-acceptable containers or conditions**

Some samples cannot be analyzed due to improper collection or due to deterioration during transfer. If any sample is deemed inadequate or needs to be rejected, you will be notified accordingly.

Samples may be rejected or deemed inadequate due to:

- Use of glass containers for frozen samples
- Containers that present leakage which affects the accuracy of the examination
- Syringes with the needle still attached
- Improper storing temperature (refrigerated, frozen, room temperature)
- Sample quantity inadequate for testing
- Wrong transportation medium
- Samples that have been collected using means whose expiration date has passed
- Improper or inadequate labeling or samples that don't match the examination requirements
- Samples with no examination request form
- Incomplete or improper examination request form
- Samples with no label
- Contaminated samples